



## Ziehm Solo FD

Versatile design meets  
latest flat-panel technology



Ziehm Solo FD. As the size of hospital and surgery center ORs decreases and equipment quantity rises, the demand for imaging systems with smaller footprints is growing. With its all-in-one design, the Ziehm Solo FD is one of the most compact C-arms for even the smallest treatment scenarios on the market – now equipped with CMOS flat-panel technology to perform a broad portfolio of applications. Versatile viewing options and new dimensions in user friendliness offer maximum flexibility in the OR to support your clinical workflow. With the latest improvements in dose regulations, the Ziehm Solo FD ensures best image quality at a minimized dose.

# 01 / Achieve significantly more details with CMOS flat-panel technology

Optimal soft tissue and bone contrast as well as high spatial resolution and a wide dynamic range are key to displaying detail-rich images of even the smallest anatomical structures. CMOS detector technology delivers on all counts, helping physicians to improve image quality.

## → CMOS flat-panel technology

Image quality and efficiency are the most important but also challenging factors in daily clinical routine. In comparison with conventional C-arms, the latest flat-panel technology CMOS achieves higher spatial resolution due to smaller pixel sizes combined with lower noise levels and a higher read-out speed at full resolution. True resolution, especially in the magnification modes, makes interpolation unnecessary. CMOS technology therefore enables improved overall efficiency.

## → Higher level of performance

The compact monoblock generator provides short, sharp pulses, producing razor-sharp images even if the patient is moving. This intelligent pulse technology also improves dose management. The flat-panel technology is unaffected by magnetic fields and enables distortion-free imaging, with no loss in image quality and more than 65,000 shades of gray.

## → Contrast-rich visualization

The Ziehm Solo FD offers a contrast-rich 19" dual display with high brightness. Even from a distance, the high-end monitors provide the physician with optimal insights by visualizing the finest details – from any angle.




Full size (7.9" x 7.9")



Magnification mode 1 (5.9" x 5.9")




Magnification mode 2 (3.9" x 3.9")



Spatial resolution phantom with more than 4.0 lp/mm visible




Conventional image intensifier



Ziehm Solo FD with flat-panel technology





Cervical spine



Stabilization of a radius fracture



Adjusting screw for stabilization of syndesmosis



Peripheral revascularization




Osteosynthesis of the clavicular

02/Ensure maximum flexibility  
with a versatile design

As space in the OR is limited, the demand for imaging systems with smaller footprints is growing. Thanks to the compact design and viewing options, the Ziehm Solo FD enables the hospital to fit every individual need.

→ Compact design

The Ziehm Solo FD is one of the smallest C-arms on the market. It comes as standard with 19" TFT dual color flatscreens mounted on an articulating monitor arm, eliminating the need for a separate monitor cart. Despite the compact design, 165 degrees of orbital movement enhances easy patient coverage and ensures a maximum of flexibility in the OR – fully counter-balanced in every position.



All-in-one design

All functions required for image capturing, processing and archiving are integrated in the C-arm, without the need for a standalone monitor cart.




Easy handling

165 degrees of orbital movement and the 34" C-arm opening ideally support your workflow.

→ Flexible configurations

In addition to the compact design of the system, three different viewing options enhance flexibility during interventions to allow the product range to suit individual needs.

These options allow you to conveniently operate the system from the Ziehm Viewing Station, the Remote Solo Center and the C-arm.



Option 1: wall- or ceiling-mounted monitors  
This space-saving configuration maximizes available space in the OR and can benefit from wireless integration.



Option 2: Remote Solo Center  
Create sufficient scope for sterile operation with the Remote Solo Center, flexibly mounted to the sides of the OR table or on a separate stand.



Ziehm Solo FD with integrated monitor  
This versatile mobile C-arm comes as standard with an integrated monitor to ensure a compact design for small ORs. Furthermore, it can be extended with three different viewing options.

Option 3: Ziehm Viewing Station  
The C-arm can be easily supplemented with an extra Viewing Station featuring a 24" landscape splitscreen monitor or dual 19" flatscreen monitors.





## 03 / Optimize process efficiency with advanced clinical workflows

In the face of time and efficiency pressure, compatible clinical workflows help to operate the C-arm in an easy and intuitive way. Unmistakable communication increases safety in the OR and optimizes efficient patient handling.

### → Wireless Freedom

WLAN allows you to store images wirelessly to the PACS from any location. With Ziehm Wireless Video, live images can be transferred to wall-or ceiling mounted monitors in real-time. Benefit from the wireless dual-plus-footswitch to increase safety due to fewer cables on the OR floor, and increase flexibility, as relevant functions such as initiating X-ray can be controlled wirelessly. Encrypted technology is used for all mentioned features.

### → Fit for the future

The Solo Center is a touchscreen with an open, modular software architecture, ensuring maximum flexibility. This interface can be easily upgraded and expanded with additional software modules without the need for hardware changes.

### → Seamless integration

The open interface, Ziehm NetPort, enables easy integration into existing IT networks. X-rays saved in DICOM 3.0 format are transferred to the PACS, and patient data can be exchanged with HIS/RIS. Data can be retrieved at any time. It can also be backed up to DVD or USB stick and printed on transparencies or paper.



Ziehm SmartEye technology mirrors the live image on the touchscreen, enabling the operator to keep track of orientation and object position.

## 04 / Work with a minimized dose with intelligent dose regulation – now in its next generation

Minimizing dose while maintaining image quality is an important goal worldwide for surgeons, their staff and patients alike. Ziehm Imaging understands and supports this by developing enhancements in different applications to help our customers reach this goal.

→ SmartDose<sup>1</sup>. Best image quality. Minimized dose.

The comprehensive concept consists of a broad, clinically proven application portfolio to answer daily challenges of low dose and high image quality. Further reduce dose significantly with dedicated SmartDose functions for pediatric surgery.<sup>2</sup> With significant dose savings, Ziehm Imaging sets the benchmark in user-friendly adjustment of dose exposure. Our latest improvements in SmartDose help to display even the smallest details of challenging anatomical areas and reduce dose with intelligent pulse regulation and optimized anatomical programs.

<sup>1</sup> In clinical practice, the use of SmartDose may reduce patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task.

<sup>2</sup> Gosch D. et al. "Influence of grid and ODDC on radiation exposure and image quality using mobile C-arms – First results", RöFo, 09/07




- LASER POSITIONING DEVICE**  
integrated in flat-panel/image intensifier and generator housing for accurate placement of C-arm without radiation
- REDUCTION OF PULSE FREQUENCY**  
pulse frequency from 1 to 30 p/s manually or automatically for lower accumulated dose
- OBJECT DETECTED DOSE CONTROL (ODDC)**  
automatically analyzes the area of interest to minimize dose and optimize image quality
- ANATOMICAL PROGRAMS – AUTOMATIC ADJUSTMENTS FOR BEST RESULTS**  
optimize dose and image quality automatically for best results in dedicated anatomical programs
- HIGH-SPEED ADR – INTELLIGENT PULSED REGULATION**  
reduce dose by taking advantage of latest innovations for faster pulse regulation
- ZAIP – ADVANCED INSIGHTS IN FAST-MOVING OBJECTS**  
benefit from latest improvement of ZAIP algorithm and filters for imaging razor-sharp details of guide-wires and even the smallest vessels
- LOW DOSE MODE**  
pediatric key for all anatomical programs to lower dose to the patient
- PREMAG**  
for exposure-free display of magnified X-ray image
- AUTOMATIC ADJUSTMENT**  
automatic adjustment for obese patients with no additional dose increase
- REMOVABLE GRID**  
for pediatric and other dose-sensitive procedures to reduce dose
- VIRTUAL COLLIMATORS**  
for exposure-free positioning of collimators



## 05 / Features at a glance

In a challenging healthcare environment, where space requirements and cost sensitivity are the greatest influences, the Ziehm Solo FD puts you forward with a wide range of versatile applications and options to fit your individual needs.

Footprint	8.6sqft	Anatomical Marking Tool (ATM)	optional
C-arm opening	34"	Ziehm Viewing Station	optional
2k x 2k CMOS technology	•	Monitors for existing ceiling support arms	optional
Touchscreen user interface	•	Ziehm NetPort	optional
Ziehm SmartEye with SmartControl	•	Interface to 2D navigation systems	optional
SmartArchive	•	Printer / DVD	optional
Color-coded scales and handles	•	Cineloop/DSA	optional
Pulsed monoblock generator	•	Wireless freedom integration (WLAN, Wireless Footswitch and Wireless Video)	optional



Color-coded handles



Wireless Footswitch and Remote Solo Center



# NATIONWIDE SERVICE COVERAGE



## 24/7

### 24 HR TOLL-FREE SERVICE AND TECHNICAL PHONE SUPPORT

- Parts inventory in the USA delivered in 24 hrs in 98 % of requests
- Comprehensive service programs from our direct service organization

Call 866-949-4346



We would be pleased to consult with you personally. Please scan to give us a call.

#### USA

Ziehm Imaging Inc.  
6280 Hazeltine National Dr.  
Orlando, FL 32822, USA  
Toll Free +1.(800)503.4952  
Phone +1.(407)615.8560  
Fax +1.(407)615.8561  
mail@ziehm.com



OrthoScan Inc.  
8212 East Evans Road  
Scottsdale AZ 85260, USA  
Toll Free +1.(866)996.0472  
Phone +1.(480)503.8010  
Fax +1.(480)503.8011  
info@orthoscan.com

#### Headquarters Germany

Ziehm Imaging GmbH  
Donaustrasse 31  
90451 Nuremberg, Germany  
Phone +49.(0)9 11.21 72-0  
Fax +49.(0)9 11.21 72-390  
info@ziehm-eu.com

#### Italy

Ziehm Imaging Srl  
Via Paolo Borsellino, 22/24  
42100 Reggio Emilia, Italy  
Phone +39.0522.61 08 94  
Fax +39.0522.61 24 77  
italy@ziehm-eu.com

#### Finland

Ziehm Imaging Oy  
Kumitehtaankatu 5  
04260 Kerava, Finland  
Phone +358.449 75 75 37  
finland@ziehm-eu.com

#### France

Ziehm Imaging S.A.R.L.  
1, Allée de Londres  
91140 Villejust, France  
Phone +33.1 69 07 16 65  
Fax +33.1 69 07 16 96  
france@ziehm-eu.com

#### Brazil

Ziehm Medical do Brasil  
Av. Roque Petroni Jr.,  
1089 cj 904  
04707-000 São Paulo, Brazil  
Phone +55.(11)3033.5999  
Fax +55.(11)3033.5997  
brazil@ziehm.com

#### China

Ziehm Medical Shanghai Co., Ltd.  
Hongqiao New Tower Centre  
Rm 06-07, 25/F  
83 Loushanguan Road  
Shanghai, P.R. China; 200336  
Phone +86.(0)21.62369903  
Fax +86.(0)21.62369916  
china@ziehm.net.cn

#### Singapore

Ziehm Imaging Singapore Pte. Ltd.  
7030 Ang Mo Kio Ave 5  
#08-53 Northstar@AMK  
Singapore 569880, Singapore  
Phone +65.639.18600  
Fax +65.639.63009  
singapore@ziehm-eu.com

