

# Suction canisters

Reusable Serres canisters have to be changed hardly ever and are extremely durable. Serres Suction canisters are certified as measuring devices with an accuracy of +/- 100 ml. Canisters are equipped with built-in brackets for mounting on walls, rail supports or trolleys. Canisters include reusable angle connectors for vacuum tubing. The compact oval-shaped 1000 ml canisters are ideal for tight places in ambulances, intensive care units and anaesthesia work stations.



## Suction canisters, transparent

Transparent canister with wide scale allow for a clear view of suctioned liquid.

Product no.	Product name	Packing size
57308	Suction canister, 1000 ml	45
57328	Suction canister, 2000 ml	24
57348	Suction canister, 3000 ml	16



## Suction canisters, blue

Blue canisters are designed for sensitive circumstances to ensure comfort of patients and families.

Product no.	Product name	Packing size
57309	Suction canister blue, 1000 ml	45
57329	Suction canister blue, 2000 ml	24



## Flushing canister

In intermittent use suction tube can be cleaned and maintained operational by suctioning small amount of water from Flushing canister.

Product no.	Product name	Packing size
57319	Flushing canister, white, 1000 ml	45



## Valves and connectors for Serres canisters

- On-off valve can be used to open and close suction. It replaces the grey angle connector.
- Serres valve is used to open and close suction when the smoke filter is used.
- T-connector is used in serial connections and can replace the Serres valve when the smoke filter is used
- The reusable Angle connector is available as a spare part if it is disposed of accidentally.



Product no.	Product name	Packing size
57958	Serres valve	10
57959	On-Off valve	10
57610	T-connector	100
57618	Angle connector	20